Spektrofotometrie

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Meas. mode	Wavelengths
Wavelength [nm]	812
Integration time [s]	0.2

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

Serial measurement

File Edit Start Reference Help Close To document

Measurement parameters
Sample series

Measurement settings Notes Table

Designation	SPECORD 40 - 2321
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Wavelengths
Wavelength [nm]	812
Integration time [s]	0.2

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Measurement settings Notes

Designation SPECORD 40 - 2321
Display Absorbance
Correction Reference
Lamp change at 320 nm
Meas. mode Wavelengths
Wavelength [nm] 812
Integration time [s] 0.2

Table

Table

Specord 40 - C:\Users\Specord40\Documents\Praktika_CZ\Kantori\X00_C.par

OK Cancel New Open Save Save as... Help Print

Settings Device Mode Accessories

Title Calibration

Operator X00

Cycle Mode
 None
 Manual
 Automatic
 Time controlled
Number 3
Interval time [s] 1

Display
 Transmittance
 Absorbance
 Energy
 Reflectance
Min: 0
Max: 3

Correction
 No
 Standard
 Reference
 Special
Open

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Serial measurement

File Edit Start! Reference Help Close To document

Measurement settings Notes Table

Designation	SPECORD 40 - 232I
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Wavelengths
Wavelength [nm]	812
Integration time [s]	0.2

Specord 40 - untitled.par

Settings Device Mode Accessories

Title: Spectrums

Operator: 00X

Cycle Mode

- None
- Manual
- Automatic
- Time controlled

Number: 3
Interval time [s]: 1

Display

- Transmittance
- Absorbance
- Energy
- Reflectance

Min: 0.0
Max: 3.0

Correction

- No
- Standard
- Reference
- Special

Open

Cycle Node: *None*

Display: *Absorbance*

Correction: *Reference*

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Title: *Spectrums*

Serial measurement

File Edit Start! Reference Help Close To document

Measurement settings Notes Table

Designation	SPECORD 40 - 2321
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Wavelengths
Wavelength [nm]	812
Integration time [s]	0.2

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Specord 40 - untitled.par

OK Cancel Save Open Save as... Help Print

Settings Device **Mode** Accessories

Meas. mode Scan Mode

Scan Mode

Range [nm] 375 - 1100

Delta lambda [nm] 1,0

Speed [nm/s] 50,0

Integration time [s] 0.02

not saved

Meas. mode: *Scan Mode*

Range [nm]: *375 -1100*

Delta lambda [nm]: *1.0*

Speed [nm/s]: *50.0*

Serial measurement

File Edit Start! Reference Help Close To document

Measurement settings Notes Table

Designation	SPECORD 40 - 2321
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Wavelengths
Wavelength [nm]	812
Integration time [s]	0.2

Klikněte na ikonu „Save as...“.

Specord 40 - untitled.par

OK Cancel New Open Save Save as... Help Print

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Serial measurement

File Edit Start! Reference Help Close To document

Measurement settings Notes Table

Designation	SPECORD 40 - 232I
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Wavelengths
Wavelength [nm]	812
Integration time [s]	0.2

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Specord 40 ... (Users\Specord40\Documents\Praktika_CZ\Kantori\00X_S.par)

OK Cancel New Open Save Save as... Help Print

Settings Device Mode Accessories

Meas. mode Scan Mode

Scan Mode

Range [nm] 375 - 1100

Delta lambda [nm] 1,0

Speed [nm/s] 50,0

Integration time [s] 0.02

saved OK

1. Klikněte na ikonu „Edit“.

Serial measurement

File Edit Start Reference Help Close To document

Measurement parameters

Sample series

Measurement settings Notes

Designation	SPECORD 40 - 2321
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Scan Mode
Range [nm]	375 - 1100
Delta lambda [nm]	1
Speed [nm/s]	50

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
Spectrums A, B and C
00X

Measurement settings

- Designation
- Display
- Correction
- Lamp change at
- Meas. mode
- Range [nm]
- Delta lambda [nm]
- Speed [nm/s]

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
Table Grafik

Designation SPECORD 40 - 2321
Display Absorbance
Correction Reference
Lamp change at 320 nm
Meas. mode Scan Mode
Range [nm] 375 - 1100
Delta lambda [nm] 1
Speed [nm/s] 50

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
Sample 1
Sample 2
Sample 3

Parameter

Designation	SPECTRUM 40-1100
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Scan Mode
Range [nm]	375 - 1100
Delta lambda [nm]	1
Speed [nm/s]	50

Transmittance Absorbance Values

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
Samples

Sample 1
Sample 2
Sample 3

Parameters Audit trail Notes

Designation SPECORD 40 - 20
Display Absorbance
Correction Reference
Lamp change at 320 nm
Meas. mode Scan Mode
Range [nm] 375 - 1100
Delta lambda [nm] 1
Speed [nm/s] 50

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




2. V podmenu vyber položku „Sample properties“.

Designat	Reference
Display	
Correction	Reference
Lamp change at	320 nm
Meas. mode	Scan Mode
Range [nm]	375 - 1100
Delta lambda [nm]	1
Speed [nm/s]	50

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

Samples

Sample 1 - XXXXX
Sample 2 - YYYYY
Sample 3 - ZZZZZ

Parameters Audit trail Notes


Designation	SPECORD 40 - 20
Display	Absorbance
Correction	Reference
Lamp change at	320 nm
Meas. mode	Scan Mode
Range [nm]	375 - 1100
Delta lambda [nm]	1
Speed [nm/s]	50

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
Samples

- Sample 1 - XXXXX
- Sample 2 - YYYYY
- Sample 3 - ZZZZZ

Parameters	Audit trail	Notes
Designation	SPECORD 40 - 20	
Display	Absorbance	
Correction	Reference	
Lamp change at	320 nm	
Meas. mode	Scan Mode	
Range [nm]	375 - 1100	
Delta lambda [nm]	1	
Speed [nm/s]	50	

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